


Reusable



ERCP-CATHETERS

SPRAY CATHETERS / WASH-OUT PROBES

PAPILLOTOMES

CYSTOSTOME / ESD-KNIVES

BASKETS

LITHOTOMY BASKETS / LITHOTRIPTOR (SET) /
BALLOON CATHETERS

STENTS

INTRODUCER SYSTEMS / STENT-SETS /
CYST DRAINAGE SETS / RING KNIFE SET /
DPOC-GUIDING CATHETER

POLYPECTOMY SNARES

BIOPSY FORCEPS

HOT-BIOPSY FORCEPS / PUNCTURE BIOPSY FORCEPS / HIGH
FREQUENCY CLAMP / MICRO-BITES

FOREIGN BODY REMOVING FORCEPS

POLYPOTOMES / POLYP TRAP NET LOOPS

BRUSHES

CYTOLOGY BRUSHES/ CLEANING BRUSHES

INJECTION NEEDLES

ASPIRATION NEEDLE

GUIDE WIRES

CUSTOMER SUPPORT INTERNATIONAL

Jutta Kasten-Heß
+49 2803 9120-118

Louiza Knops
+49 2803 9120-119

Djamila Knops
+49 2803 9120-120

ADMINISTRATION

Ulrike Ackermann
+49 2803 9120-117

Elke Großbongardt
+49 2803 9120-116

QUALITY MANAGER

Robert Koch
+49 2803 9120-122

COLOUR CODES AND SYMBOLS



— Disposable products



— Reusable products



— Single-lumen products



— Double-lumen products



— Triple-lumen products



— Shape-stabilizer



0.035 — Suitable for guide wire 0.035 inch



0.025 — Suitable for guide wire 0.025 inch



0.018 — Suitable for guide wire 0.018 inch



250 — Suitable for double- and single-balloon enteroscopy, 250 cm



280 — Suitable for double- and single-balloon enteroscopy, 280 cm



120 — Suitable for bronchoscopy, 120 cm


We reserve the right to make technical changes and to deviate from the depicted images.

ERCP-CATHETERS



ERCP-CATHETERS

partly with inside metal ring, radiopaque, Luer-Lock connection, mandrin removable




Article no.	Catheter tip	Diameter mm	Length cm	For guide wire inch
0120100	standard	1,6	215	0.025
0120200	filiform	1,6	215	0.025
0130200	filiform	1,8	215	0.035
0120210	filiform / with inside metal ring	1,6	215	0.025

Figure shows article REF 0120100, 0130200

ERCP-CATHETERS

metal tip, Luer-Lock connection, mandrin removable




Article no.	Catheter tip	Diameter mm	Length cm	For guide wire inch
0120410	filiform metal tip	1,6	215	0.025
0120510	bottle-shaped metal tip	1,6	215	0.025
0120610	conical metal tip	1,6	215	0.025

Figure shows article REF 0120410, 0120510, 0120610

ERCP-CATHETERS

tapered, partly with inside metal ring, radiopaque,
guide wire lockable, integrated mandrin




Article no.	Catheter tip	Diameter mm	Length cm	For guide wire inch
0121110	standard / with inside metal ring	2,3 > 1,6	215	0.025
0130101	standard	2,3 > 1,8	215	0.035
0131110	standard / with inside metal ring	2,3 > 1,8	215	0.035
0120201	filiform	2,3 > 1,6	215	0.025
0120211	filiform / with inside metal ring	2,3 > 1,6	215	0.025
0130201	filiform	2,3 > 1,8	215	0.035
0130211	filiform / with inside metal ring	2,3 > 1,8	215	0.035

Figure shows article REF 0120201, 0130101

ERCP-CATHETERS

tapered, with metal tip, guide wire lockable, integrated mandrin



Article no.	Catheter tip	Diameter mm	Length cm	For guide wire inch
0120311	round metal tip	2,3 > 1,6	215	0.025
0130311	round metal tip	2,3 > 1,8	215	0.035
0120411	filiform metal tip	2,3 > 1,6	215	0.025
0130411	filiform metal tip	2,3 > 1,8	215	0.035
0120511	bottle-shaped metal tip	2,3 > 1,6	215	0.025
0130511	bottle-shaped metal tip	2,3 > 1,8	215	0.035
0120611	conical metal tip	2,3 > 1,6	215	0.025
0130611	conical metal tip	2,3 > 1,8	215	0.035


Figure shows article REF 0130311, 0130411, 0130511, 0130611

SPRAY CATHETERS AND WASH-OUT PROBES



SPRAY CATHETERS

mandrin removable



mandrin removable

Article no.	Diameter mm	Length cm	
110101210	1,8	120	
0101310	2,3	215	

Figure shows article REF 0101310

WASH-OUT PROBE



Article no.	Catheter tip	Diameter mm	Length cm
110120100	without metal head, sheath with ring marking	1,6	120

PAPILLOTOMES CYSTOSTOME ESD-KNIVES

COLOUR CODES AND SYMBOLS



— Disposable products



— Reusable products



— Single-lumen products



— Double-lumen products



— Triple-lumen products



— Shape-stabilizer



0.035 — Suitable for guide wire 0.035 inch



0.025 — Suitable for guide wire 0.025 inch



0.018 — Suitable for guide wire 0.018 inch



250 — Suitable for double- and single-balloon enteroscopy, 250 cm




280 — Suitable for double- and single-balloon enteroscopy, 280 cm



120 — Suitable for bronchoscopy, 120 cm


We reserve the right to make technical changes and to deviate from the depicted images.

PAPILLOTOMES



PAPILLOTOMES

standard for guide wire




Article no.	Special features	Feather strengthened	For guide wire inch	Cutting length mm	Cutting direction	Diameter mm	Length cm
023031	5 mm guiding tip	yes	0.035	30	11 h	2,3 > 1,8	215
023034	5 mm guiding tip	yes	0.035	20	11 h	2,3 > 1,8	215

Figure shows article REF 023031

PAPILLOTOMES

filiform for guide wire




Article no.	Special features	Feather strengthened	For guide wire inch	Cutting length mm	Cutting direction	Diameter mm	Length cm
023051	7 mm radiopaque tip	yes	0.035	30	11 h	2,3 > 1,8	215
023651	7 mm radiopaque tip	yes	0.035	20	11 h	2,3 > 1,8	215

Figure shows article REF 023051

PAPILLOTOMES

precut for guide wire




Article no.	Feather strengthened	For guide wire inch	Cutting length mm	Cutting direction	Diameter mm	Length cm
023011	yes	0.035	20	11 h	2,3 > 1,8	215
023021	yes	0.035	30	11 h	2,3 > 1,8	215

Figure shows article REF 023011

PAPILLOTOME, COMBINED

feather strengthened, needle retractable



Article no.	Special features	Needle length mm	Feather strengthened	Cutting length mm	Cutting direction	Diameter mm	Length cm
020101	needle retractable	10	yes	30	11 h	1,6	215

PAPILLOTOMES


needle
knife



Article no.	Papillotome type	Cutting length mm	Diameter mm	Length cm
020121	needle	6	1,6	215
020131	knife	6	1,6	215

Figure shows article REF 020121


HF-NEEDLE



Article no.	Papillotome type	Needle length mm	For guide wire inch	Diameter mm	Length cm
023121	needle	6	0.035	2,3 > 1,8	215

PAPILLOTOME

for B-II-stomach



Article no.	For guide wire inch	Maximum cutting arc mm	Diameter mm	Length cm
023161	0.035	20	2,3	215

COLOUR CODES AND SYMBOLS



— Disposable products



— Reusable products



— Single-lumen products



— Double-lumen products



— Triple-lumen products



— Shape-stabilizer



— Suitable for guide wire 0.035 inch



— Suitable for guide wire 0.025 inch



— Suitable for guide wire 0.018 inch



— Suitable for double- and single-balloon enteroscopy, 250 cm



— Suitable for double- and single-balloon enteroscopy, 280 cm



— Suitable for bronchoscopy, 120 cm

We reserve the right to make technical changes and to deviate from the depicted images.

LITHOTOMY BASKETS LITHOTRIPTOR (SET) BALLOON CATHETERS


03

BASKETS



LITHOTOMY BASKET

6 wires, with nozzle for contrast medium, suitable for lithotripsy




injection port for contrast medium

Article no.	Basket expanded dimension mm		Basket type	Basket wire	Diameter mm	Length cm
	Opening width	Opening length				
03123704	35	50	slightly turned	twisted	2,6	215

LITHOTOMY BASKET


4 wires, with nozzle for contrast medium, **suitable for lithotripsy**



Article no.	Basket expanded dimension mm		Basket type	Basket wire	For guide wire inch	Diameter mm	Length cm
	Opening width	Opening length					
03144434	20	30	round	twisted	0.035	2,6	215

LITHOTOMY BASKET

4 to 8 wires with nozzle for contrast medium, **not** suitable for lithotripsy




injection port for contrast medium

Article no.	Basket expanded dimension mm		Basket wire	Diameter mm	Length cm
	Opening width	Opening length			
03113608	30	45	twisted	2,6	215

LITHOTOMY BASKETS

4 wires, twisted with nozzle for contrast medium, suitable for lithotripsy



Article no.	Basket expanded dimension mm		Basket type	Basket wire	For guide wire inch	Diameter mm	Length cm
	Opening width	Opening length					
03112324	15	25	diamond shaped	twisted	0.025	1,8	215
03112424	20	30	diamond shaped	twisted	0.025	1,8	215
03112524	25	40	diamond shaped	twisted	0.025	1,8	215
03112624	30	45	diamond shaped	twisted	0.025	1,8	215
03112724	40	55	diamond shaped	twisted	0.025	1,8	215
03114434	20	30	diamond shaped	twisted	0.035	2,6	215
03114534	25	40	diamond shaped	twisted	0.035	2,6	215
03114634	30	45	diamond shaped	twisted	0.035	2,6	215
03114734	40	55	diamond shaped	twisted	0.035	2,6	215
03114834	45	60	diamond shaped	twisted	0.035	2,6	215

Figure shows article REF 03114634

COLOUR CODES AND SYMBOLS



— Disposable products



— Reusable products



— Single-lumen products



— Double-lumen products



— Triple-lumen products



— Shape-stabilizer



— Suitable for guide wire 0.035 inch



— Suitable for guide wire 0.025 inch



— Suitable for guide wire 0.018 inch



— Suitable for double- and single-balloon enteroscopy, 250 cm



— Suitable for double- and single-balloon enteroscopy, 280 cm




— Suitable for bronchoscopy, 120 cm

We reserve the right to make technical changes and to deviate from the depicted images.

POLYPECTOMY SNARES




POLYPECTOMY SNARE, SYMMETRICAL



Article no.	Opening size width x height mm	Wire type	Isolated tip	Diameter mm	Length cm
052011202	10 x 20	twisted	no	2,3	230
052011222	15 x 30	twisted	no	2,3	230
052011212	20 x 40	twisted	no	2,3	230
052011232	30 x 55	twisted	no	2,3	230
052031202	10 x 20	twisted	yes	2,3	230
052031212	20 x 30	twisted	yes	2,3	230
052031232	30 x 55	twisted	yes	2,3	230

Figure shows article REF 052031212


POLYPECTOMY SNARE, SYMMETRICAL



Article no.	Opening size width x height mm	Wire type	Diameter mm	Length cm
052011225	20 x 40	monofil	2,3	230
052011235	25 x 50	monofil	2,3	230

Figure shows article REF 052011225


POLYPECTOMY SNARE, SYMMETRICAL WITH FLAT FEATHER



Article no.	Opening size width x height mm	Wire type	Isolated tip	Diameter mm	Length cm
052011206	10 x 24	twisted with flat feather	no	2,3	230
052011226	20 x 40	twisted with flat feather	no	2,3	230
052011236	30 x 60	twisted with flat feather	no	2,3	230
052031226	20 x 40	twisted with flat feather	yes	2,3	230
052031236	30 x 60	twisted with flat feather	yes	2,3	230

Figure shows article REF 052031226

POLYPECTOMY SNARE, HEXAGONAL



Article no.	Opening size width x height mm	Wire type	Isolated tip	Diameter mm	Length cm
052011234	30 x 50	twisted	no	2,3	230
052031234	30 x 50	twisted	yes	2,3	230

Figure shows article REF 052031234


POLYPECTOMY SNARE, ASYMMETRICAL



Article no.	Opening size mm	Isolated tip	Diameter mm	Length cm
051011121	25	no	1,8	160
052011211	15	no	2,3	230
052011221	25	no	2,3	230
052031211	15	yes	2,3	230
052031221	25	yes	2,3	230

Figure shows article REF 052031221

POLYPECTOMY SNARE, ASYMMETRICAL



Article no.	Opening size mm	Isolated tip	Diameter mm	Length cm
052010248	40	no	2,3	230
052010268	60	no	2,3	230
052030248	40	yes	2,3	230
052030268	60	yes	2,3	230

Figure shows article REF 052030248


POLYPECTOMY SNARE, SHARK TEETH, SYMMETRICAL



Article no.	Opening size width x height mm	Wire type	Diameter mm	Length cm
052051225	20 x 40	monofil, with 4 teeth	2,3	230
052061235	25 x 50	monofil, with 4 teeth	2,3	230

Figure shows article REF 052051225

POLYPECTOMY SNARE "VARIO-SNARE", SYMMETRICAL



Article no.	Opening size width x height mm	Wire type	Isolated tip	Diameter mm	Length cm
052091232	approx. 15 x 33 (step 1) or approx. 27 x 50 (step 2)	twisted	yes	2,3	230

COLOUR CODES AND SYMBOLS



— Disposable products



— Reusable products



— Single-lumen products



— Double-lumen products



— Triple-lumen products



— Shape-stabilizer



0.035 — Suitable for guide wire 0.035 inch



0.025 — Suitable for guide wire 0.025 inch



0.018 — Suitable for guide wire 0.018 inch



250 — Suitable for double- and single-balloon enteroscopy, 250 cm



280 — Suitable for double- and single-balloon enteroscopy, 280 cm



120 — Suitable for bronchoscopy, 120 cm

We reserve the right to make technical changes and to deviate from the depicted images.

BIOPSY FORCEPS
HOT-BIOPSY FORCEPS
PUNCTURE BIOPSY FORCEPS
HIGH FREQUENCY CLAMP
MICRO-BITES

06

BIOPSY FORCEPS



BIOPSY FORCEPS

fenestrated, oval spoonshaped mouth, tapered



Article no.	Diameter mm	Length cm
1106140210	1,8	120
061402101	1,8	160
061402104	1,8	230
1106150210	2,2	120
061512101	2,2	160
061512105	2,2	180
061512104	2,2	230
061602101	2,6	160
061602104	2,6	230

Figure shows article REF 061402104

BIOPSY FORCEPS

fenestrated, oval spoonshaped mouth, with spike, tapered




Article no.	Diameter mm	Length cm
1106140211	1,8	120
061402114	1,8	230
061512111	2,2	160
061512115	2,2	180
061512114	2,2	230
061602111	2,6	160
061602114	2,6	230

Figure shows article REF 061402114

BIOPSY FORCEPS

oval spoons shaped mouth, toothed



Article no.	Diameter mm	Length cm
061003004	0,75	230
061103003	1,0	230

Figure shows article REF 061103003

HOT-BIOPSY FORCEPS

oval spoonshaped mouth




Article no.	Diameter mm	Length cm
062502004	2,2	230
062602004	2,6	230

Figure shows article REF 062502004

The hot-biopsy forceps are isolated with a PTFE coat and do have a lateral injection port.

HIGH FREQUENCY CLAMP

model Dr. Straumann



Article no.	Shape of high frequency clamp	Diameter mm	Length cm
062600004	2 grasping arms with cross grooves	2,2	230

The high frequency clamp is isolated with a PTFE coat and do has a lateral injection port.

BIOPSY FORCEPS

controllable, oval spoonshaped mouth, fenestrated



controllable

injection port

Article no.	Diameter mm	Length cm
1106460210	2,6	120
064602105	2,6	180

Figure shows article REF 064602105

COLOUR CODES AND SYMBOLS



— Disposable products



— Reusable products



— Single-lumen products



— Double-lumen products



— Triple-lumen products



— Shape-stabilizer



0.035 — Suitable for guide wire 0.035 inch



0.025 — Suitable for guide wire 0.025 inch



0.018 — Suitable for guide wire 0.018 inch



250 — Suitable for double- and single-balloon enteroscopy, 250 cm



280 — Suitable for double- and single-balloon enteroscopy, 280 cm




120 — Suitable for bronchoscopy, 120 cm

We reserve the right to make technical changes and to deviate from the depicted images.

FOREIGN BODY REMOVING 07
FORCEPS
POLYPOTOMES
POLYP TRAP NET LOOPS

FOREIGN BODY REMOVING FORCEPS



FOREIGN BODY REMOVING FORCEPS


PTFE tube, with lateral injection port



Article no.	Opening dimension mm	Number of arms	Diameter mm	Length cm
0720123	30	3	2,5	230

FOREIGN BODY REMOVING FORCEPS

metal tube



Article no.	Opening dimension mm	Number of arms	Diameter mm	Length cm
1107101	20	3	1,8	120
130710120	20	3	1,8	250

Figure shows article REF 1107101

FOREIGN BODY REMOVING FORCEPS

metal tube



Article no.	Opening dimension mm	Forceps shape	Diameter mm	Length cm
1107107	7,8	oval spoon shaped, toothed	1,8	120
1107207	7,2	oval spoon shaped, toothed	2,2	120

Figure shows article REF 1107107

FOREIGN BODY REMOVING FORCEPS

metal tube



Article no.	Opening dimension mm	Forceps shape	Diameter mm	Length cm
0721020	8	long alligator mouth, flat jaws with cross and lengthways grooves	2,2	230

MTW

FOREIGN BODY REMOVING FORCEPS “PEG-FORCEPS”

metal tube



Article no.	Opening dimension mm	Diameter mm	Length cm
0721230	8	2,6	180

FOREIGN BODY REMOVING FORCEPS, TISSUE FIXATION FORCEPS

metal tube



Article no.	Opening dimension mm	Diameter mm	Length cm
0721240	8	2,6	185

FOREIGN BODY REMOVING FORCEPS, GRASPING FORCEPS FOR NUTRITION PROBES


metal tube



Article no.	Opening dimension mm	Diameter mm	Length cm
0731243	10	2,5	320

POLYPOTOMES

PTFE tube, with lateral injection port




Article no.	Opening dimension mm	Number of arms	Diameter mm	Length cm
0720223	30	3	2,5	230
0734423	30	4	2,5	230
0735523	30	5	2,5	230

Figure shows article REF 0720223

POLYPOTOME, STATIONARY

PTFE tube, with lateral injection port



Article no.	Opening dimension mm	Number of arms	Diameter mm	Length cm
0720224	25	3	2,6	230

COLOUR CODES AND SYMBOLS



— Disposable products



— Reusable products



— Single-lumen products



— Double-lumen products



— Triple-lumen products



— Shape-stabilizer



0.035 — Suitable for guide wire 0.035 inch



0.025 — Suitable for guide wire 0.025 inch



0.018 — Suitable for guide wire 0.018 inch



250 — Suitable for double- and single-balloon enteroscopy, 250 cm



280 — Suitable for double- and single-balloon enteroscopy, 280 cm



120 — Suitable for bronchoscopy, 120 cm

We reserve the right to make technical changes and to deviate from the depicted images.

BRUSHES
CYTOLOGY BRUSHES
CLEANING BRUSHES

BRUSHES



CYTOLOGY BRUSH, WITH METAL BALL



Article no.	Brush diameter mm	Diameter ball mm	Diameter mm	Length cm
130812201	2,0	1,8	1,8	250

CLEANING BRUSHES, DOUBLE-SIDED



Article no.	Brush diameter mm	Working channel mm	Diameter mm	Length cm
0803102	3,0	2,0	1,8	160
0803202	3,0	2,0	1,8	230
0805102	5,0	2,8	1,8	160
0805202	5,0	2,8	1,8	230

Figure shows article REF 0805102

CLEANING BRUSHES, ONE-SIDED



Article no.	Brush diameter mm	Working channel mm	Diameter mm	Length cm
0803101	3,0	2,0	1,8	160
0803201	3,0	2,0	1,8	230
0805101	5,0	2,8	1,8	160
0805201	5,0	2,8	1,8	230

Figure shows article REF 0805101

COLOUR CODES AND SYMBOLS



— Disposable products



— Reusable products



— Single-lumen products



— Double-lumen products



— Triple-lumen products



— Shape-stabilizer



— Suitable for guide wire 0.035 inch



— Suitable for guide wire 0.025 inch



— Suitable for guide wire 0.018 inch



— Suitable for double- and single-balloon enteroscopy, 250 cm




— Suitable for double- and single-balloon enteroscopy, 280 cm




— Suitable for bronchoscopy, 120 cm

We reserve the right to make technical changes and to deviate from the depicted images.

INJECTION NEEDLES



INJECTION NEEDLES



Article no.	Cannula diameter mm	Extended cannula length mm	Diameter mm	Length cm
0910518218	0,5	5	1,8	180
0910718218	0,7	5	1,8	180
0910723223	0,7	5	2,3	230

Figure shows article REF 09 10 51 82 18

INJECTION NEEDLES

with lateral nozzle to wash-out



Article no.	Cannula diameter mm	Extended cannula length mm	Diameter mm	Length cm
0910718523	0,7	5	2,3	180
0910723523	0,7	5	2,3	230

Figure shows article REF 0910718523

INJECTION NEEDLE

with metal tube



Article no.	Cannula diameter mm	Extended cannula length mm	Diameter mm	Length cm
0910723123	0,7	5	2,3	230

COLOUR CODES AND SYMBOLS



— Disposable products



— Reusable products



— Single-lumen products



— Double-lumen products



— Triple-lumen products



— Shape-stabilizer



— Suitable for guide wire 0.035 inch



— Suitable for guide wire 0.025 inch



— Suitable for guide wire 0.018 inch



— Suitable for double- and single-balloon enteroscopy, 250 cm



— Suitable for double- and single-balloon enteroscopy, 280 cm



— Suitable for double- and single-balloon enteroscopy, 550 cm




— Suitable for bronchoscopy, 120 cm



— Extra long


We reserve the right to make technical changes and to deviate from the depicted images.

GUIDE WIRES



GUIDE WIRES, NITINOL

kink-resistant, 50 mm flexible platinum tip




Article no.	Guidewire tip	Color version coiled	Diameter inch	Length cm
001800	straight platinum tip	blue	0.018	400
002200	straight platinum tip	green-black	0.025	400
002255	straight platinum tip	green-black	0.025	550
002300	straight platinum tip	green-black	0.025	460
003200	straight platinum tip	yellow-black	0.035	400
003226	straight platinum tip	yellow-black	0.035	260
003500	angled platinum tip, "J-Tip"	yellow-black	0.035	400
003446	straight platinum tip, "Soft-Tip"	yellow-blue	0.035	460
003600	angled platinum tip, "J-Soft-Tip"	yellow-blue	0.035	400
003615	angled platinum tip, "J-Soft-Tip"	yellow-blue	0.035	150

Figure shows article REF 003200


GUIDE WIRE “PLATIN-STAR”

kink-resistant, Nitinol, with scale



Article no.	Guide wire tip	Color version coiled	Diameter inch	Length cm
003800	straight platinum tip	yellow-black	0.035	460

This guide wire is equipped with 25 visible endoscopic and X-ray platinum rings and 7 visible endoscopic ring marks



MTW[®]

Endoskopie Manufaktur



MTW-Endoskopie W. Haag KG | Goldbergstraße 18 | D-46487 Wesel

Tel +49 2803 9120-0 | Fax +49 2803 9999788 | info@mtw-endoskopie.com | www.mtw-endoskopie.com
